

## INTRODUCTION

The ScriptCap™ 2'-O-Methyltransferase Kit prepares cap 1-RNA from any source of cap 0-RNA. In part, cap 1 methylation serves, to increase the translation efficiency of the mRNA.<sup>1</sup> ScriptCap™ 2'-O-Methyltransferase transfers a methyl group from the donor molecule S-adenosyl-methionine (SAM) to the 2'-O position of the penultimate nucleotide of a cap 0 RNA (m<sup>7</sup>G[5']ppp[5']NpN...) to synthesize RNA with a cap 1 structure (m<sup>7</sup>Gppp[m<sup>2'-O</sup>]NpNpN...). The cap 0 RNA can be produced by enzymatically capping uncapped RNA using the ScriptCap™ m<sup>7</sup>G Capping System or by *in vitro* transcription of a DNA template in the presence of a dinucleotide cap analog (e.g., m<sub>2</sub><sup>7,3'-O</sup>GpppG; e.g., using a MessageMAX™ T7 ARCA-Capped Message Transcription Kit). Cap 1 RNA can also be synthesized from uncapped RNA in a single reaction mixture that contains both the ScriptCap m<sup>7</sup>G Capping System and ScriptCap 2'-O-Methyltransferase plus SAM.

One ScriptCap 2'-O-Methyltransferase Kit reaction methylates 60 µg of 5'-Cap 0 capped RNA.

## MATERIALS

### Materials Supplied



Do not store in a frost-free freezer. Do not store at -70°C.

| ScriptCap™ 2'-O-Methyltransferase Kit Contents (25 reactions)  |                |
|--|----------------|
| Kit Component  | Reagent Volume |
| ScriptCap™ 2'-O-Methyltransferase, (blue cap)<br>100 U/µl, in 50% glycerol, 50 mM Tris-HCl, pH 7.5, 100 mM NaCl,<br>1 mM dithiothreitol (DTT), 0.1 mM EDTA and 0.1% Triton® X-100. | 100 µl         |
| 10X ScriptCap™ Capping Buffer<br>0.5 M Tris-HCl, pH 8.0, 60 mM KCl and 12.5 mM MgCl <sub>2</sub>   | 250 µl         |
| 20 mM S-adenosyl-methionine (SAM)  | 65 µl          |
| ScriptGuard™ RNase Inhibitor, 40 U/µl<br>in 50% glycerol, 50 mM Tris-HCl, pH 7.5, 100 mM NaCl, 10 mM<br>DTT, 0.1 mM EDTA and 0.1% Triton X-100.                                    | 65 µl          |
| RNase-Free Water   | 2 x 1.4 ml     |

Inquire about custom kit sizes at 608-442-6484 or sales@cellscript.com.

### Materials Required, but not Supplied

- Cap 0 RNA substrate
- Materials or kits for purification of the RNA product

## SPECIFICATIONS

### Unit Definition

One unit of ScriptCap 2'-O-Methyltransferase methylates one picomole of a control Cap 0 RNA in 1 hour at 37°C under standard assay conditions.

### Contaminating Activity Assays

All components of the ScriptCap 2'-O-Methyltransferase Kit are free of detectable RNase and DNase activities.

## BEFORE YOU START: IMPORTANT TIPS FOR OPTIMAL CAPPING

### ◆ SAM:

SAM slowly degrades over time at room temperature and above. Keep thawed SAM solutions on ice at all times.

### ◆ Cap 0 RNA Source:

RNA should be purified and resuspended in RNase-Free Water. **Do not resuspend the RNA in an EDTA-containing solution.**

- a) Cap 0-RNA produced using the capping enzyme-based ScriptCap m<sup>7</sup>G Capping System (sold separately): Directly add the ScriptCap 2'-O-Methyltransferase to the ScriptCap m<sup>7</sup>G Capping System reaction either simultaneously or sequentially without prior reaction clean-up.
- b) Cap 0-RNA generated using a dinucleotide cap analog in an *in vitro* transcription reaction: Cleanup the RNA prior to treatment with ScriptCap 2'-O-Methyltransferase. Purify the RNA by your preferred method. The method chosen should remove residual proteins and unincorporated nucleotides from the RNA.

### ◆ RNA Secondary Structure:

Some RNA transcripts can form stable secondary structures (homodimers and hairpins) involving the 5'-most nucleotides of the transcript severely limiting access of the 5'-penultimate nucleotide to the ScriptCap 2'-O-Methyltransferase. In order to increase the enzymatic efficiency on such RNAs, use longer or hotter heat denaturation conditions than those noted in the protocol. Times and temperatures required will vary.

### ◆ Poly(A)-Tails:


If the cap 1-RNA requires subsequent 3'-poly(A)-tailing, using CELLSCRIPT's A-Plus Poly(A) Polymerase (sold separately) allows the user to skip RNA purification prior to poly(A)-tailing (see the A-Plus Poly(A) Polymerase product literature for details). Capped and tailed RNA must be purified prior to use in RNA transfection experiments.


## PROCEDURE

## A. Synthesis of Cap 1-RNA from Cap 0-RNA

1. The protocol below was designed for use with 50-60 µg of cap 0-RNA. Combine the following reagents:


| Standard ScriptCap 2'-O-Methyltransferase Kit Reaction (step 1) |        |
|---|--------|
| Component   | Amount |
| RNase-Free Water  | x µl   |
| Cap 0-RNA, 50-60 µg   | ≤81 µl |
| Total Volume  | 81 µl  |


 Heat-denaturation of the RNA is an optional step, but it is strongly recommended for RNAs which have not previously been characterized for their ease of enzymatic capping.

 Only heat-denature the RNA and water components. **Do not** include any other reagent in this step.

2. Incubate at 65°C for 5-10 minutes, then transfer to ice.
3. While the heat-denatured RNA is cooling on ice, prepare a "Cocktail" of the following reaction components together in a separate tube.


| Standard ScriptCap 2'-O-Methyltransferase Kit Reaction (step 3) |        |
|---|--------|
| Component   | Amount |
| 10X ScriptCap Capping Buffer                                    | 10 µl  |
| 20 mM SAM   | 2.5 µl |
| ScriptGuard RNase Inhibitor                                     | 2.5 µl |
| ScriptCap 2'-O-Methyltransferase, 400 Units                     | 4 µl   |
| Total Volume  | 19 µl  |

 A white precipitate may form in the 10X ScriptCap Capping Buffer upon storage. To dissolve it, heat the tube at 37°C for 5 minutes and mix thoroughly.

 Keep the thawed stock of SAM on ice.

4. Combine the Cap 0 RNA Solution from Step 1 with the Cocktail from Step 3.

| Standard ScriptCap 2'-O-Methyltransferase Kit Reaction (step 4) |        |
|---|--------|
| Component   | Amount |
| Heat-denatured cap 0-RNA (from step 1)                          | 81 µl  |
| Cocktailed reaction components (from step 3)                    | 19 µl  |
| Total Reaction Volume   | 100 µl |

 The efficiency of 2'-O-methylation can be lower if the RNA 5' end is structured. If your RNA is not completely 2'-O-methylated, we recommend increasing the incubation time to 2 hours. Also, since the concentration of methylation sites for a given mass is higher for small RNA than for large RNA, increase the reaction time for small RNA. For example, we suggest to increase the reaction time to 2 hours if your RNA is <730 nucleotides in length.


5. Incubate at 37°C for 30 minutes.
6. The Cap 1 RNA can now be purified, or it can be 3' polyadenylated without purification by adding the reaction mixture directly to an A-Plus™ Poly(A) Polymerase reaction (sold separately). Purification of poly(A)-tailed Cap 1 RNA is recommended prior to use for RNA transfection.


**B. Simultaneous Capping and 2'-O-Methylation to Synthesize Cap 1 RNA from Uncapped RNA**

1. You need to purchase the ScriptCap m<sup>7</sup>G Capping System in addition to ScriptCap 2'-O-Methyltransferase in order to synthesize cap 1 RNA from uncapped RNA using the following protocol. This protocol was designed for use with 50-60 µg of uncapped RNA.

Combine the following reagents:


| Simultaneous Capping and 2'-O-Methylation (step 1) |        |
|--|--------|
| Component  | Amount |
| RNase-Free Water                                   | x µl   |
| <i>In vitro</i> transcribed uncapped RNA, 50-60 µg | ≤67 µl |
| Total Volume                                       | 67 µl  |


 Heat-denaturation of the RNA is an optional step, but it is strongly recommended for RNAs which have not previously been characterized for their ease of enzymatic capping.


 Only heat-denature the RNA and water components. **Do not** include any other reagent in this step.

2. Incubate at 65°C for 5-10 minutes, then transfer to ice.
3. While the heat-denatured RNA is cooling on ice, prepare a "Cocktail" of the following reaction components together in a separate tube.

| Simultaneous Capping and 2'-O-Methylation (step 3) |        |
|--|--------|
| Component  | Amount |
| 10X ScriptCap Capping Buffer                       | 10 µl  |
| 10 mM GTP *  | 10 µl  |
| 20 mM SAM  | 2.5 µl |
| ScriptGuard RNase Inhibitor                        | 2.5 µl |
| ScriptCap 2'-O-Methyltransferase (100 U/µl)        | 4 µl   |
| Total Volume                                       | 29 µl  |

 Do not include the ScriptCap Capping Enzyme in this mix.


 Keep the thawed stock and diluted SAM solutions on ice.

 A white precipitate may form in the 10X ScriptCap Capping Buffer upon storage. To dissolve it, heat the tube at 37°C for 5 minutes and mix thoroughly.

4. **Just prior to starting the reaction**, add the ScriptCap Capping Enzyme to the Cocktail from Step 3 and then combine this with the uncapped RNA solution from Step 1.

| Simultaneous Capping and 2'-O-Methylation (step 4) |        |
|--|--------|
| Component  | Amount |
| Cocktailed reaction components (from step 3)       | 29 µl  |
| ScriptCap Capping Enzyme (10 U/µl)*                | 4 µl   |
| Heat-denatured RNA (from step 1)                   | 67 µl  |
| Total Reaction Volume                              | 100 µl |

\* The 10 mM GTP in Step 3 and ScriptCap Capping Enzyme in Step 4 are components of the ScriptCap m<sup>7</sup>G Capping System.

 The efficiency of 2'-O-methylation can be lower if the RNA 5' end is structured. If your RNA is not completely 2'-O-methylated, we recommend increasing the incubation time to 2 hours. Also, since the concentration of methylation sites for a given mass is higher for small RNA than for large RNA, increase the reaction time for small RNA. For example, we suggest to increase the reaction time to 2 hours if your RNA is <730 nucleotides in length.

5. Incubate at 37°C for 30 minutes.
6. The Cap 1 RNA can now be purified, or it can be 3' polyadenylated without purification by adding the reaction mixture directly to an A-Plus™ Poly(A) Polymerase reaction (sold separately). Purification of poly(A)-tailed Cap 1 RNA is recommended prior to use for RNA transfection.

## TROUBLESHOOTING

| Symptom   | Solution   |
|---|--|
| <b>Low capping or 2'-O-methylation efficiency</b> | RNA to be treated with ScriptCap 2'-O-Methyltransferase should be purified and resuspended in RNase-free water. Do not resuspend the RNA in an EDTA-containing solution.   |
|   | Prior to 2'-O-methylation, purify the input RNA using a method that removes residual proteins, contaminants and unincorporated nucleotides and/or cap analogs.   |
|   | Verify that ScriptGuard RNase Inhibitor was added to the reaction.   |
|   | SAM slowly degrades at room temperature and above. Keep SAM solutions on ice at all times.   |
|   | Increase the reaction incubation time. For example, up to 3 hours at 37°C.   |
| <b>White precipitate in reaction buffer</b>       | Some RNAs form stable structures (e.g., homodimers, hairpins) at the 5' end, limiting access by Capping Enzyme or 2'-O-Methyltransferase. Analyze the sequence and increase the RNA denaturation temperature to above the $T_m$ (e.g., to 65°C for 20 min, 75°C for 10 min, 85°C for 5 min). If the 5' end is highly structured, it might be necessary to modify the 5' end sequence using molecular biology techniques. Often this can be accomplished by making a single point mutation within the first 5 bases of the DNA template for the RNA transcript (non-coding region). |
|   | Incubate the reaction buffer at 37°C for 5 minutes then mix thoroughly to dissolve the precipitate.  |
|   | Do not store the kit at -70°C.   |

## RELATED PRODUCTS

- A-Plus™ Poly(A) Polymerase Tailing Kit
- INCOGNITO™ SP6  $\Psi$ -RNA Transcription Kit
- INCOGNITO™ T7 5mC- &  $\Psi$ -RNA Transcription Kit
- INCOGNITO™ T7  $\Psi$ -RNA Transcription Kit
- INCOGNITO™ T7 ARCA 5mC- &  $\Psi$ -RNA Transcription Kit
- MessageMAX™ T7 ARCA-Capped Message Transcription Kit
- ScriptCap™ m<sup>7</sup>G Capping System
- ScriptGuard™ RNase Inhibitor
- SP6-Scribe™ Standard RNA IVT Kit
- T7 mScript™ Standard mRNA Production System
- T7-FlashScribe™ Transcription Kit
- T7-Scribe™ Standard RNA IVT Kit

## REFERENCE

1. Kuge, H. et al., (1998) Nucl. Acids Res. 26, 3208.

The performance of this product is guaranteed for one year from the date of purchase.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

A-Plus, AmpliCap, AmpliCap-Max, CELLSCRIPT, INCOGNITO, MessageMAX, mScript, ScriptCap, ScriptGuard, SP6-Scribe, T7-FlashScribe, T7-Scribe and 2-Way are trademarks of CELLSCRIPT, Madison, Wisconsin.

The Purchaser of this product agrees to the TERMS AND CONDITIONS posted on CELLSCRIPT's website: <http://www.cellscript.com>.

© 2017 CELLSCRIPT, All rights reserved.